



## ZIUZ IRIS

The IRIS checks 13 different parameters such as shape, size and colour thus greatly reducing DAA packing errors – it has a very low error rate of 1.33 sachets per million.

The IRIS is very compact and takes up the space of a printer and thus is ideal for pharmacies with space restraints.

The ZiuZ IRIS takes a photo of each dosage bag for data comparison. These photos are saved in a central storage, so that they can be easily retrieved if necessary.

The IRIS uses advanced visual intelligence technology and various parameters; it analyses the contents of each medication sachet, identifying errors on screen which can then be checked and corrected by the Pharmacist.